








Upper Anteriors | Antérieures maxillaires | Anteriores superiores | Anteriori superiori






Ovoid
Ovale
Forma de pala
Ovale





Triangular
Triangulaire
Triangular
Triangolare




Rectangular
Rectangulaire
Rectangular
Rettangolare



Square
Carrée
Cuadrado
Quadrato



Lower Anteriors | Antérieures mandibulaires | Anteriores inferiores | Anteriori inferiori



* Long (O47 or L37 with longer tooth neck) F: Longue (O47 et L37 avec collet plus long) ES: Largo (O47 o L37 con cuello más largo) I: Lungo (O47 e L37 con colletto allungato)

VITA MFT Anterior Antérieures Anteriores Anteriori	UJ / No. mm	O40 40.0	O44 43.9	O49 49.2	T41 40.5	T43 43.3	T46 46.1	T49 48.8	R42 42.1	R45 44.5	R47 46.6	S47 46.8	S50 49.7
	LJ / No. mm	L33 32.6	L34 34.0	L40 39.5	L33 32.6	L34 34.0	L37 36.8	L40 39.5	L34 34.0	L37 36.8	L37 36.8	L40 39.5	L40 39.5
VITA MFT Posterior Postérieures Posteriores Posteriori	UJ / No. mm	PU29 28.8	PU29 28.8	PU33 33.0	PU29 28.8	PU31 30.6	PU31 30.6	PU33 33.0	PU29 28.8	PU31 30.6	PU33 33.0	PU33 33.0	PU33 33.0
	LJ / No. mm	PL29 29.6	PL29 29.6	PL33 33.6	PL29 29.6	PL31 31.1	PL31 31.1	PL33 33.6	PL29 29.6	PL31 31.1	PL33 33.6	PL33 33.6	PL33 33.6

VITA MFT®


shades | couleurs | colores | colori:

16 x VITA SYSTEM 3D-MASTER® (0M1, 0M2, 0M3, 1M1, 2L1.5, 2M1, 2M2, 3M1, 3L1.5, 3L2.5, 3M2, 3R2.5, 3M3, 4L1.5, 4M2, 5M1)


16 x VITA classical A1-D4® (A1, A2, A3, A3.5, A4, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4)

Posterior | Postérieures | Posteriores | Posteriori

PU29 – PL29



PU31 – PL31



PU33 – PL33

